() Riester ri-thermo® N

Clinical tympanic thermometer

Is the new ergonomic infrared thermometer for uncomplicated use in medical practices, hospitals and in the home. Its innovative infrared technology not only allows **ri-thermo® N** to measure the body temperature of children and adults easily in a matter of seconds but also means that it can measure the temperature of surfaces, rooms and bathwater. The temperature values are displayed in a range from 0 to 100 °C or on the corresponding scale in ° Fahrenheit. Its large LCD display with backlighting makes the values easy to read. Even in the dark. Displays up to 12 stored temperatures. The **ri-thermo® N** sounds a handy acoustic signal when it has finished measuring and if it detects a fever, has an energy-saving automatic switch-off function and comes with a hygienic storage container with probe covers. **ri-thermo® N** is of the usual high **Riester** quality, is clinically validated and is supplied with batteries.



| Measurement range: | 0 °C to 100.0 °C (32.0 °F to 212.0 °F) |
|---------------------------------|---|
| Measuring accuracy: | ± 0.2 °C, 32.0 - 42.2 °C (± 0.4 °F, 89.6 - 108.0 °F) |
| Display: | Liquid crystal display with 0.1 °C (0.1 °F) display units |
| Acoustic: | a. The device is switched on (ON) and ready to measure: 1 short beep. b. The measuring process has finished: 1 long beep. c. System error or fault: 3 short beeps |
| Memory: | a. The last measured value appears when the device is switched on.b. 12 measurements can be retrieved in memory mode. |
| Operating temperature: | 5 °C - 40 °C (41.0 °F - 104 °F) |
| Storage/ transport temperature: | -25 °C - +55 °C (-13 °F - 131 °F) |
| Automatic switch-off: | About 1 minute after the last measurement. To prevent misuse, the device switches itself off automatically when the START button is pressed and held down for more than 15 sec. |
| Battery: | CR2032 BATTERY (X1) – at least 1000 measurements. |
| Size: | 153 mm (L) x 31 mm (W) x 40 mm (H) |
| Weight: | 53 g (with battery), 50 g (without battery) |
| Standards: | Meets ASTM E-1965 requirements |
| | |
| No. 1805 | ri-thermo® N with 25 probe covers |
| No. 12700 | Disposable probe covers, pack of 25 pcs. |
| No. 12701 | Disposable probe covers, pack of 100 pcs. |
| No. 12702 | Disposable probe covers, pack of 1000 pcs. |

Please visit riester.de for more information

